# EXELYA® 愛尔雅

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## INDICATIONS FOR USE

The device is a super-oxidized cleanser intended to help fight acne.

Through fluid flowing through the lesion, it mechanically removes biologic and inert materials, reduces the amount of microorganisms and contributes to a moist environment, so supporting the creation of an ideal microenvironment, enabling the body to perform its physiological healing process.

## HANDLING

Non special handling precautions are required.

## STORAGE

The device should be stored in its original sealed container between 5°C and 25°C. It is non-flammable. Special storage precautions are not needed.

## STABILITY

Until opened, the stability of the device is assured up to expiration date (see bottle for year and month of expiration). Upon opening, note the date on the bottle in the section "Date Opened". Stability of the unused solution is assured for up to 30 days after opening.

## DIRECTIONS FOR USE

EXELYA® is ready for use with no mixing or dilution required. Apply EXELYA® directly on the affected area:

- Open the bottle by unscrewing the cap.
- Open the plastic sachet containing the spraying cap by cutting the upper side.
- Extract the spraying cap from the sachet, by taking it out from the cap, paying attention not to touch the straw.
- Insert the straw into the bottle.
- Screw the spraying cap on the bottle neck.
- Before using the product for the first time, spray on air to start the pump, till you see a nebulized solution.
- Apply EXELYA® directly on the injured zone, by spraying it, in order to profusely cover the affected area, let it dry, repeat the application.
- Apply the treatment twice a day: in the morning and in the evening.
- Do not unscrew the spraying cap after its installation; do not apply again the screw cap to the bottle after the installation of the spraying cap.
- Close with the overcap just after use.
- Use the device within 30 days after the first opening.

For the treatment administration, if the patient is going to use make up, he/she should do so only after the product is dried for at least two minutes after its application. If the patient has to shave, he/she should do so before the application of the product.

## WARNINGS AND POINTS TO NOTE

Not for injection or consumption.

Keep out of the reach and sight of children.

Keep away from food and beverages.

Do not ingest, inhale or put the device in contact with eyes.

Do not dilute the product before use.

Use the device within 30 days after the first opening.

Do not unscrew the spraying cap after its installation; do not apply again the screw cap to the bottle after the installation of the spraying cap.

Close with the overcap just after use.

Store the medical device protected from light and heat, between 5°C and 25°C, in its carton box, in a dry place.

Do not use the device if the packaging is damaged or if the security sealer is already away or broken.

Do not use the device after the expiry date.

Do not use in case of known or suspected individual hypersensitivity to any component of the product (Hypochlorous acid (HOCl), Cl<sub>2</sub> and Hypochlorite ion (ClO<sup>-</sup>)).

The use in pregnancy and during breastfeeding has not been tested.

Local irritation or itching may occur; with the prolonged use sensitization may happen. In case, stop treatment and ask for medical advice.

In case of adverse events, stop the treatment and inform your physician.

In case of severe acne, or if you have any treatment concerns, consult a health care professional.

If acne does not show signs of healing, consult health care professional.

Always use a fresh solution, do not re-use the solution which have already been in contact with wounds. The re-use can cause lack of efficacy, infection, sensitization/irritation phenomena.

The use of the product in combination with other devices/drugs/cosmetics has not been tested.

Do not put the bottle entrance or the spray cap in contact with the patient body and/or dressing in order to avoid cross-contamination.

Dispose the solution according to the local dispositions for expired medicinal/medical device products.

## COMPOSITION

Super-oxidized water solution. It contains hypochlorous acid (HClO).

Use the device within 30 days after the first opening.

Store the medical device protected from light and heat, between 5°C and 25°C, in its carton box, in a dry place.

**Date of issue of the latest revision of the instruction for use:**  
December 2016

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